

TEKNİK ÖZELLİKLER

İçerik : Mikro gözenekli non-woven rayon fiber kumaş
Yapışma Gücü : En az 120 gf/cm
Su Buharı Geçirgenliği : En az 500 g/m²
TS 3957 e uygundur
Raf Ömrü : 2 yıldır

TECHNICAL PROPERTIES

Contain : The micro pored non-woven rayon fabric
Adhesive Power : Min 120 gf/cm
Water Steam Permeability : Min. 500 g/m²
It is proper to TS 3957
Shelf Lifetime : 2 years

EBAT / SIZES	KOLİ İÇİ(Adet) / BOX IN (pcs)	KOLİ EBATLARI (cm) / BOX SIZES (cm)
5 x 1,25	200	32 x 21 x 30
5 x 2,5	120	34 x 31 x 20
5 x 5	100	37 x 24 x 32
10 x 2,5	100	34 x 31 x 20